

ALMANET

January 2015

Magnesium - Powder MGP 100/300µm

Application:	Pyrotechnic								
Composition:	Made from Mg 99,9 % min.								
Granulation:	Mg powder 100-300 µm Typical sieve analysis to JEL 200 sieve inlet: DIN ISO 3310-1 <table><tr><td>> 300 µm</td><td>0%</td></tr><tr><td>> 200 µm</td><td>30 – 60%</td></tr><tr><td>> 100 µm</td><td>90 – 100%</td></tr><tr><td>< 56 µm</td><td>0%</td></tr></table>	> 300 µm	0%	> 200 µm	30 – 60%	> 100 µm	90 – 100%	< 56 µm	0%
> 300 µm	0%								
> 200 µm	30 – 60%								
> 100 µm	90 – 100%								
< 56 µm	0%								
Apparent Density:	0,75 – 0,85 g / cm ³								
Special Properties:	Due to a special production procedure a low oxide content is achieved								
Packing:	UN-certified 210-l steel drums								

